

### Eyewash bottles oculav NIT®

SÖHNGEN®

#### Execution

Ready-to-use sterile solution for instant rinsing of eyes, the conjunctival sac and affected skin where the accident occurred. The buffer solution binds and neutralises, especially after acid or alkaline burns. Bottle can be kept and used for three years provided unopened. Maintenance-free. Tested in accordance with DIN EN 15154-4.



### Technical Data

Product number	1093009 195
DIN	15154-4
Manufacturer Part Number	2010080
EAN Code	4260725760535
Commodity Code	30049000



Hommel Hercules Werkzeughandel GmbH & Co. KG

Heidelberger Str. 52 | 68519 Viernheim | +49 (0) 6204 739-0 | info@hommel-hercules.com