

Instant eyewash oculav NIT® box

Execution

Ready-to-use sterile solution for instant rinsing of eyes, the conjunctival sac and affected skin where the accident occurred. The buffer solution binds and neutralises, especially after acid or alkaline burns. Bottle can be kept and used for three years provided unopened. Maintenance-free. Tested in accordance with DIN EN 15154-4.

Delivery

oculav-NIT® box, 4 x 250 ml bottle oculav NIT® sterile solution in orange plastic box with wall bracket, 260 x 160 x 80 mm, also with immediate use set for eye injuries with DuOcul eye compresses, eye patch, SÖHNGEN® silk plaster strips, DERMOTEKT® compress and disposable protective gloves

SÖHNGEN®



Technical Data

Product number	1093009 196
DIN	15154-4
Colour	Orange
Manufacturer Part Number	2010004
EAN Code	4250108801504
Commodity Code	30065000



Hommel Hercules Werkzeughandel GmbH & Co. KG

Heidelberger Str. 52 | 68519 Viernheim | +49 (0) 6204 739-0 | info@hommel-hercules.com