

# PRODUCT DATA SHEET

**Eyewash bottles** Sodium chloride solution

### Execution

Eyewash bottle with sterile 0.9% sodium chloride solution for rinsing out foreign objects such as dust, dirt and splinters in the event of mechanical injury. The quick and easy-to-use bottles are suitable for mobile use on site, e.g. in tool cases or first aid cases, and for refilling wall stations art. no. 1093009 186 and art. no. 1093009 188. Flushing times: 200-ml bottle approx. two min., 500-ml bottle approx. five min. Can be kept and used for three years provided unopened. Maintenance-free. Tested in accordance with DIN EN 15154-4.

#### Notes

Flushing times: 200-ml bottle approx. two min., 500-ml bottle approx. five min.





## **Technical Data**

Product number	1093009 183
DIN	15154-4
Manufacturer Part Number	4691
EAN Code	5701684469125
Commodity Code	33079000

