

### Eye emergency station

#### Execution

Eye emergency station in dust-protected wall-mounted box for rapid first aid in cases of eye-related accidents. Contents: 200 ml eyewash bottle with sterile 4.9% phosphate buffer solution for neutralisation of chemical substances, 500 ml eyewash bottle with sterile 0.9% sodium chloride solution for rinsing out foreign objects such as dust, dirt and splinters in case of mechanical injury. A eye-related emergency plan is also included. Flushing times: 200 ml bottle approx. two min., 500 ml bottle approx. five min. Bottles can be kept and remain usable for three years provided unopened. Maintenance-free. Tested in accordance with DIN EN 15154-4.

#### Delivery

Wall box including 1 x 200 ml eyewash bottle with sterile 4.9% phosphate buffer solution and 1 x 500 ml eyewash bottle with sterile 0.9% sodium chloride solution

#### Notes

Flushing times: 200 ml bottle approx. 2 min., 500 ml bottle approx. 5 min.

plum



### Technical Data

Product number	1093009 188
DIN	15154-4
Colour	Green
Manufacturer Part Number	4789
EAN Code	5701684478912
Commodity Code	33079000



Hommel Hercules Werkzeughandel GmbH & Co. KG

Heidelberger Str. 52 | 68519 Viernheim | +49 (0) 6204 739-0 | info@hommel-hercules.com